

HEARTSCREEN 112C-1

3/6/12-channel diagnostic ECG



Main features:

- 12-lead simultaneous measurement, monitoring
- 3/6/12-channel, 112mm wide thermal printing
- Foil keys and alphanumeric keyboard on touchscreen display for patient data input
- Last 20 patient data storage
- 8-hours operating from a battery with a 2-hours charging time
- Built-in power supply
- 150 recordings storage capacity
- Diagnosis -, rhythm analysis - and rhythm event printing program
- Automatic and manual lead selection
- Programmable recording modes
- Automatic ECG curve positioning
- Economic printing mode
- Quick paper load
- USB or WiFi connection to PC
- Connection to hospital information system
- Exercise ECG option on PC
- Cardiological data management software on PC



MEDICON SERVICES

Destination to Medical Solution

416/1, Ground Floor, Mokki Masjid Road, West Rampura, Dhaka-1219

Mobile: 01711823264

Email: mediconservices416@gmail.com

Phone: +88 02-55128410

cmohsinkamal@yahoo.com

www.mediconservices.com



Made in Hungary