


Human

Made in Germany

Diagnostics Worldwide

Biochemistry Analyzer

HumaLyzer 3000



Model Name

Humalyzer 3000

Spectrophotometer type

Filter Photometer

Optical Configuration

Single beam with continuously rotating filter wheel

Monochromatic or bichromatic reading

8 filter positions

Usable Spectral Range

330 to 670 nm

System Procedures

Open and by stored menu

Calculating Modes

Absorbance

Single standard

Differential samples

Factor Mode

Differential samples

Multi Standard Mode (up to 7 standards)

Multi Standard % abs(up to 7 standards)

Kinetic Mode (Consecutively or simultaneously Batch)

By Factor or By Standard



MEDICON SERVICES

Destination to Medical Solution

416/1, Ground Floor, Mokki Masjid Road, West Rampura, Dhaha-1219

Mobile: 01711823264

Email: mediconservices416@gmail.com

Phone: +88 02-55128410

cmohsinkamal@yahoo.com

www.mediconservices.com

