

# MIKASA HF100H


Origin Japan



## HF100H

High power general-purpose x-ray unit.  
Ideal choice for nursing homes.

The HF100H is compact and lightweight enabling easy transportation. It is capable of taking exposures of most body parts using maximum 100 kV and maximum 40mA high frequency technology. The unit has 5 memory stations for storing x-ray exposure techniques. It is compatible with all major DR digital systems. It is recommended for hospitals, mobile clinics, military field clinics and nursing homes.



Input voltage	: 220-230V, AC single phase, 50/60Hz
Output	: 2kW
kV range	: 40 – 100kVDC
mA range	: 20 – 40mA
Timer	: 0.02 sec. – 4.0 sec.
Focal spot value	: 1.2
Dimension	: 286mm (W) x 400mm (L) x 230mm (H)
Weight	: 16.8kg
Others	: rotating collimator, 5 x memory of output settings.



**MEDICON SERVICES**

416/1 Ground Floor, Mokki Mosjid Road,  
West Rampura, Dhaka-1219.

Phone: +8802-55128410  
mediconservices416@gmail.com

Mobile: +8801711823264.  
www.mediconservices.com