

Point 200 HD



Economic with High Quality



For every dentist

- The Point 200HD is a user friendly device with clear panoramic image.

Point 200 HD®

Only Panoramic



• Excellent rate X-ray absorption

.....
CCD CMOS absorbs 35% of X-ray, so it is needed more X-ray. However Point 200 HD's Smart Sensor absorbs 100% of X-ray, and it is shown high quality image with lower X-ray dosage.


• High Definition image

.....
Provides more clear image by smart sensor which is applying new technology.

• More Easy to do a Patient Positioning

.....
Patient positioning is assisted by cross positioner system which accurately indicates the correct anatomical positioning points.

● Adoption Smart Sensor



Easily shown clear and high quality image when the luminous body absorbed X-ray, the each cell converts Electronic signal without any spread and overlap.




The Paniramic image in previous time



The Paniramic image currently

Point 200 HD	
Type	CMOS
Sensor Type	Smart Sensor(CdTe)
Pixel Size	48 um
Total Active Area	3,125 x 100 mm
A/D Convert	12 bit
Scan Time	17 sec

Dimension



📍 416/1, Ground Floor, Mokki Masjid Road, West Rampura, Dhaka-1219, Bangladesh.

☎ +880 1711-823264 ✉ mediconservices416@gmail.com

☎ +88 02-55128410 📧 cmohsinkamal@yahoo.com

📘 mediconservices 🌐 www.mediconservices.com