



PX-500FC

**High Class X Ray Imaging
Made Affordable**

PROFESSIONAL IMAGING INC.



Specification	PX-500FC
Output Power (As per IS:7620)	40KW
RadioGraphy Parameters	500mA-100-kV kV Range 45-100 kV (varying type) Small Focus (50, 100,200) Large Focus (300,400, 500)
Exposure Time	0.2 to 5 Sec.
Time Steps	24 Steps 32 Steps Through C 24 Steps 32 Steps Through C
Exposure Release	From Control console with Two step Switch or from examination unit
Mains Voltage Compensation	Provided (354-447 volts)
Filament Stabilisation	Provided
Over Load Protection	Provided, Blocking for any hazardous combination
Tube Head	Canon Tube Model: E7239X Dual Focus 20/40 kW
Collimeter	Manual with electric timer
HV Cable	Pair of HV cable
Power Supply 3 Phase	440V+10%-15%, 50Hz 63A,
Main Resistance	<_ .36 ohm



MEDICON SERVICES

Destination to Medical Solution

📍 5A, Level 5, Ahsanullah Anima Vista, 30 Topkhana Road, Dhaka-1000.

☎ +880 1711-823264

✉ mediconservices416@gmail.com

📘 mediconservices

☎ +88 02-55128410

cmohsinkamal@yahoo.com

🌐 www.mediconservices.com