

Wall Mounted Dental X-ray Machine

Model: JYF-10B



Specification

Name	Dental X Ray Machine New JYF-10B Wall Mounted Suction Unit
Power supply	AC220V±10%,50HZ, 1KVA
Tube voltage	70KV
Tube current	8mA
Focus size	0.8mm
Total filtration	2.5mmAL
Exposure time	0.2-4s
Leak radiation	Outside one meter≤0.002Mg/h(National standard:0.25Mg/h)
Package size	93×36×48(cm)
Net weight	22kg
Gross weight	30kg